

## 16Second\_Anatomical\_Stitch Study

GE Locus Ultra (“vct”)

- **Login** uid:sysadmin, pwd:guest
- Choose “Scan” **Mode** from LHS main screen menu
- Use **Study**: 16Second Anatomical Stitch 86mm Table Translation
- Choose The **Protocol** you want to use (Must be selected before creating exam in the next step) (For now, choices are “16 second Anatomical” and “16 second 80kV”)
- [FIX] Cradle choices
- **Create Exam** Named: YYYYMMDD MyObject; click “Create” button
- **Landmark the object:**
  - o Enter scanner room, mount pink, foam, custom designed bed (for 5-segment scans) on scanner table [FIX - DESCRIBE AND TAKE PICTURE OF POSITIONING]
  - o Secure object on custom-designed bed using tape or velcro
  - o Turn on alignment lasers - press “Toggle Lasers” soft key button on STC monitor (on the front face of the scanner)
  - o Align internal laser marker with custom bed markers (between bed-edge slots outlined in black marker – as shown FIX TAKE PICTURE) using “Fast Advance”, “Advance Table”, “Retract Table”, and “Fast Retract” soft kwy buttons on STC monitor
  - o Set the “zero” scan position - press “Set” on the internal (or external) Landmarks soft key button on the STC monitor
- **Acquire Brightfield And Darkfield Images**
  - o highlight desired scan sequence from the list of scan Names on main scanner console
  - o click “Execute” (bottom right corner of main menu)
    - “Scan” field in lower right corner must read “Bright Dark”
  - o press “Start Scan” button on main console keyboard after it glows green
  - o Wait for scan to complete (approx. 40 seconds)
- **Acquire Segment 1 Scan**
  - o Highlight desired scan sequence
    - eg. free 16spr 100vpr 102mm 680rows, 80kVp, 50mA
  - o In the sequences table, ensure that the:
    - “Start” field shows “0.0” and
    - “Bright/Dark” field reads “Completed”
  - o click “Execute” (bottom right corner of main menu)
    - “Scan” field in lower right corner must read “Sequence”
  - o press “Start Scan” button on main console keyboard after it glows green
  - o Wait for scan to complete (approx. 90 seconds)
  - o NB: Make note of Scan Series Number for all acquisitions

- **Acquire Segment 2 Scan**

- Highlight desired scan sequence (in this case, the same one)
  - eg. *free 16spr 100vpr 102mm 680rows, 80kVp, 50mA*
- In the sequences table, change the “Start” field to “-86.0”.
- In the sequences table, ensure that the
  - “Start” field shows “-86.0” and
  - “Bright/Dark” says “Completed”
- click “*Execute*” (bottom right corner of main menu)
  - “Scan” field in lower right corner must read “Sequence”
- press “*Start Scan*” button on main console keyboard after it glows green
- Wait for scan to complete (approx. 90 seconds)
- NB: Make note of Scan Series Number for all acquisitions
- Repeat for Segment 3 (or additional) Scans or Click “*Close Exam*” button

- **Acquire Segment 3 (or additional) Scan(s)**

- Highlight desired scan sequence (in this case, the same one)
  - eg. *free 16spr 100vpr 102mm 680rows, 80kVp, 50mA*
- In the sequences table, change the “Start” field to “-172.0”.

*Uber-Geek Table Translation Notes:*

- *Segment 1 translation = 0.0 mm*  
*Segment 2 translation = -86.0 mm*  
*Segment 3 translation = -172.1 mm*  
*Segment 4 translation = -258.1 mm*  
*Segment 5 translation = -344.1 mm*  
*Segment 6 translation = -516.1 mm*  
*Segment 7 translation = -602.2 mm (likely not possible)*
  - *Useable axial FOV per acquisition volume (559 voxels)*  
*= 559 voxels x 0.153884 mm/voxel*  
*= 86.02156 mm*
  - *Table displacement for N<sup>th</sup> scan is (N-1) x 86.02156 mm*
    - *Round to the nearest integer mm*  
*Scanner does not support no sub-mm bed translation*
  - *Maximum landmark-able table travel is 531.8 mm*
- 
- In the sequences table, ensure that the
    - “Start” field shows “-172.0” and
    - “Bright/Dark” says “Completed”
  - click “*Execute*” (bottom right corner of main menu)
    - “Scan” field in lower right corner must read “Sequence”
  - press “*Start Scan*” button on main console keyboard after it glows green
  - Wait for scan to complete (approx. 90 seconds)
  - Click “*Close Exam*” button

## Transfer Data for Offline Reconstruction and Post-Processing

- Open a terminal on a UNIX/LINUX computer (likely one of: darkness, hampton or mingus)
  - Check for available space on a data drive (such as /darkness/data, /hampton/data, or /hampton/data2) using UNIX 'df' command
    - `df -k` or `df -h`
  - 
  - Make appropriate directory tree under target directory (/darkness/data/vct, /hampton/data
  - Use secure copy 'scp' to copy all
    - Eg. `scp -r /vol/data/vct/4392[1-2]`  
[youruserID@hampton.imaging.robarts.ca:/hampton/data/vct/bio\\_specimen/morphology/moth/mothra/20110819/](mailto:youruserID@hampton.imaging.robarts.ca:/hampton/data/vct/bio_specimen/morphology/moth/mothra/20110819/)
  - Delete data off of vct scanner only after the reconstruction and post-processing are complete
- **Offline Reconstruction & Volume Stitching**
- Run the recon and stitching script (both are integrated into the script) from the parent directory to where the scan series directories reside:
    - `cd /hampton/data/vct/bio_specimen/morphology/moth/mothra/20110819/` (you should see the acquired scan series number directories within this parent directory)
    - `/mingus/data2/spollman/bin/proclog.vCT_Recon_Stitch_16sec.sh`
      - bed translation distance is editable (TABLE\_TRANSLATION\_MM=86)
      - reconstructed voxel spacing is editable (VOXEL\_SPACING=0.153884)

## POSSIBLY ADD TO SOP LATER

- choose protocols to add to study:
  - 16 second Anatomical
  - 16 second 80kv
- choose specimen:
  - Stitch\_Object