

GE Locus Ultra ("VCT") uCT Scanner

Operating Notes SOP

notes:

P.I. = David Holdsworth x24154

uCT Facility Manager = Joseph Umoh x24180 (jumoh@imaging.robarts.ca)

Troubleshooting Team: Hristo Nikolov, Chris Norley, Steve Pollmann (x24059, x24027)

items in this document are either buttons to push or icons to click or **items to select** from a menu

Start Scanner

Troubleshooting pre-start check options:

- Confirm that recon cluster has booted-up successfully - look on cold room Recon cluster monitor: "beowolf login:"

- Confirm that data drive (6TB ProRAID - /vol/data/vct) in cold room is powered on and running

- Confirm that STC Reset toggle switch is set to Enable, unless you want to prevent gantry motion for troubleshooting

- **NOTE: SET STC Reset TO "DISABLE" TO PREVENT GANTRY MOTION,**

IF EXTERNAL MONITOR OR OTHER EXTERNAL CABLES ARE CONNECTED TO DAS COMPUTER

- Confirm that STC has turned on - look for front panel LED reflection from your hand if window side cover is off for troubleshooting

- Confirm that DAS has turned on - click on "creepy-guy" icon in bottom-right corner of Fedora menu on VCT console

- View System error log - click on

Double-click Locus Ultra Console icon on desktop

Note: can also type **VCTConsole** from a LINUX Terminal to see errors when booting STC

user name: sysadmin

passwd: guest

choose **scan application**

(then wait until system comes up - will take several minutes until 11 icons show up in left column of console)

scan icon

Setup Experiment Protocols

Study: **16Second_Anatomical_Stitch_86mm_Table_Translation**
Specimen: **Stitch_Object**
Orientation: (prone = belly down; supine = back down)
Protocol: **16 second 80kV**

Create Exam: **2011mar18_LargeWaterbathScan**

Scan: **Release Interlock** Execute

Scan: **Scanner Warmup Sequence** Execute

Start Scan (large green hardware button on keyboard - will be pressed many times during warmup!)

Acquire Brightfield and Darkfield Scans - for each sequence

Sequences Name: **free_16spr_1000vpr_102mm_680rows 80kv 50mA** (repeat for each sequence in protocol that will be used)

Scan: **Bright Dark** Execute

Start Scan (after the button is lit)

Sequences Name: **free_16spr_1000vpr_102mm_680rows 120kv 20mA** (repeat for each sequence in protocol that will be used)

Scan: **Bright Dark** Execute

Start Scan (after the button is lit)

Sequences Name: **free_16spr_1000vpr_102mm_680rows 140kv 13mA** (repeat for each sequence in protocol that will be used)

Scan: **Bright Dark** Execute

Start Scan (after the button is lit)

Perform CT Acquisitions

Sequences Name **free_16spr_1000vpr_102mm_680rows 120kv 20mA**

Run Go Positional

Start Scan (after the button is lit)

wait for scan to complete and image to appear in Right monitor -

Sequences Name **free_16spr_1000vpr_102mm_680rows 120kv 20mA**
(drag scan extent box to desired position - superimposed on image)

End Positional

Sequences Name **free_16spr_1000vpr_102mm_680rows 120kv 20mA**

Scan: **Sequence**Execute

Start Scan (after the button is lit)

(repeat for scans at other energies or other specimens)

Perform CT Reconstructions

Troubleshooting Guide

20111110

Console computer would not boot (known problem - Solution: keep rebooting)

- power down system
- power off main breaker